



# i-STAT hs-TnI CONVERSION GUIDE

This guide provides assistance to customers transitioning from an existing testing method to testing with *i-STAT hs-TnI*, including components for conversion, testing options and a list of helpful resources. Supply orders, pricing, and service agreements can be finalized with an Abbott *i-STAT* sales specialist.

## INTENDED USE FOR i-STAT hs-TnI CARTRIDGE

The *i-STAT hs-TnI* cartridge with the *i-STAT 1* System is intended for use in the in vitro quantification of cardiac troponin I (cTnI) in whole blood or plasma samples in point of care or clinical laboratory settings. The *i-STAT hs-TnI* cartridge with the *i-STAT 1* System is intended to be used as an aid in the diagnosis of myocardial infarction (MI).



## TRANSITION CHECKLIST

### HARDWARE

- *i-STAT 1* analyzer
- *i-STAT* printer (if used)
- Downloader (downloader cable connections)

### CHECKS

- Analyzer software is updated
- Printer paper is available (if needed)
- Downloaders are programmed
- Downloaders are connected to network

**NOTE:** Before starting, check all equipment is available and operational. Contact *i-STAT* Technical Support if reprogramming of downloaders is needed. Connectivity may not be applicable for all facilities.

### CONSUMABLES

- Cartridges
- Controls & Cal Ver
- Transfer device
- Sample collection device

### CHECKS

- Storage for cartridges, controls, and Cal Ver according to IFU
- Identify the most appropriate transfer devices
- Check dates on collection devices (do not use expired collection devices)
- Verify collection devices are according to IFU

**NOTE:** New product codes/part numbers should be loaded into ordering system to eliminate ordering delays.

## COMPLIANCE

### METHOD EVALUATION & THE i-STAT hs-TnI CARTRIDGE

Follow internal policy and procedures approved by laboratory director; direct follow-up questions to accrediting organization.

**Abbott Method Evaluation support** | The method evaluation is not a manufacturer's requirement and specific details or information related to the above activities may be obtained from your accreditation or regulatory organizations.

**i-STAT Supply Calculator** | *i-STAT* supplies will be calculated according to the number of analyzers and method evaluation plan. Optional quality control plan, training cartridges and any additional cartridges for interface testing will be calculated at the customer's request.

**Data Collection Worksheets** | Worksheets are available in Abbott's online [i-STAT hs-TnI Resource Center](#). If a method evaluation report is needed, follow the instruction provided in the data collection worksheets. Worksheet requirements must be met for an Abbott representative to generate a report.

## PROCEDURES

Review procedures and identify changes needed. Next, obtain required signatures to revised procedures.

### PROCEDURE UPDATES TO CONSIDER

- Ranges: reportable, sex-specific, reference intervals
- Review necessary changes to elevated, actionable, and critical ranges
- Modify quality control plan procedures in accordance with IQCP (if applicable)
- Update reagent and control material storage conditions
- Visit the [Instructions for Use](#) page in the *i-STAT System Support* area of the [globalpointofcare.abbott](http://globalpointofcare.abbott) website.

### INDIVIDUALIZED QUALITY CONTROL PLAN (IQCP)

- Re-evaluate risk assessment based on the addition of a new test to the plan
- Ensure controls are done in accordance with IQCP
- Re-evaluate QC Plan and Quality Assessment in accordance with changes to risk assessment
- Lab director review and approval is required for all changes

Additional information about updates /changes are available from the [i-STAT hs-TnI Resource Center](#).

### PROFICIENCY TESTING (PT)

Laboratories performing non-waived testing may enroll in an approved PT program for regulated analytes listed in Subpart I Proficiency, Testing Programs for Nonwaived Testing. Proficiency Testing is used to evaluate a laboratory's performance.

**NOTE:** Consult with your regulatory agency to determine if PT is required for hs-TnI as an unregulated analyte. Abbott has identified API and WSLH providers as most compatible with *i-STAT hs-TnI* assay. Contact your proficiency test provider to determine available material compatible with *i-STAT hs-TnI*.

#### RESOURCES:

- Subpart I of the CLIA regulations: [Proficiency Testing programs for nonwaived testing](#)
- CMS Brochure: [Proficiency Testing and PT Referral brochure 10.15.24](#)
- API: Cardiac Markers 5 (program 140/920) [American Proficiency Institute 2025 Catalog](#)
- WSLH: Cardiac Markers (Item PT01260) [Wisconsin State Laboratory of Hygiene \(WSLH\) 2025 Catalog](#)

**NOTE:** Catalogs listed above are subject to change. Check proficiency test provider's website for most recent catalog releases.

## CONNECTIVITY

### NETWORKING & INTERFACING (applies if connectivity is available or planned. Otherwise, skip this section)

Alert IT, facilities management, and biomed resources regarding the implementation of a new analyte:

- Validate device connectivity
- Update Data Manager, LIS and HIS with new test information
  - Changes can include: **Test code** **Test name** **Test Ranges** **Action ranges**
- Verify that order sets, power plans, and test batteries are updated for *i-STAT hs-TnI* testing.

### DE CUSTOMIZATION CHANGES

Contact [i-STAT Technical Support](#) or your data management vendor for options to update your DE Customization Workspace for *i-STAT hs-TnI* compatibility.

- Links to [DE customization videos](#) and customization instructions (found in the *i-STAT 1 System Manual*) are linked from the [i-STAT hs-TnI Resource Center](#) at [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott): visit **Support > product-installation-training** to navigate to the *hs-TnI Resource Center* web page.

**NOTE:** Consult with your data management vendor to update to the latest DE version. If testing is not currently interfaced and will not be interfaced in the conversion, the interfacing section does not apply.

## EDUCATION & TRAINING

### CLINICAL AWARENESS

Inform clinicians and department leads of changes, which may include:

- Test name
- Ordering location
- Reportable ranges
- Clinical differences between conventional cTnI and high-sensitivity TnI assays (refer to *i-STAT cTnI vs hs-TnI Product Differences Flashcard*)
- Clinical support communication

### TRAINING

Resources available for your staff when converting to *i-STAT hs-TnI* testing:

- Conversion education from *i-STAT* support team
- Real-time, interactive online trainings session: Download the training schedule from the [Online Training Schedules](#) page (found at [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott) > *i-STAT Customers* menu > *i-STAT System Support* page > [Training & Education](#) page)

## VISIT THE *i-STAT* hs-TnI RESOURCE CENTER TO GET STARTED

Abbott has created a dedicated [online resource center](#) with materials you need to make the transition, including:

- Conversion Guide (this document)
- Product Differences Flashcard
- Guidance On Saving hs-TnI Samples *i-STAT hs-TnI* IFU
- Interactive Remote Training schedules
- *i-STAT Learning System* (hs-TnI cartridge)
- Schedule for *Meet the Experts* sessions
- Data Collection Worksheet (must be used per the instructions in order for Abbott to generate a report)
- Suggestions for procedure updates
- Suggestions for IQCP updates

**NOTE:** Not all documents are required for transition, but are available for use according to your facility's needs.



▶ SCAN THE CODE TO VISIT THE *i-STAT* hs-TnI RESOURCE CENTER

## NEED HELP? CONTACT OUR SUPPORT TEAM AT:



### MEET THE EXPERTS:

Experienced & expert staff to assist with questions & provide guidance  
Calendar available [HERE](#)



### VIA MONITORED EMAIL:

[hs-TnISupport@abbott.com](mailto:hs-TnISupport@abbott.com)

DEDICATED TO *hs-TnI* INQUIRIES WITH RESPONSE-GOAL WITHIN 48 BUSINESS HRS.



### VIA PHONE:

1-844-256-9531

DEDICATED TO *i-STAT* *hs-TnI* INQUIRIES ONLY

For in vitro diagnostic use only. For complete product information, visit [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott)

This material is for use in the U.S. only.

© 2024 Abbott All rights reserved | *i-STAT* is a trademark of Abbott

Abbott | Point of Care Diagnostics | 400 College Road East | Princeton, NJ 08540 | (609) 454-9000

*i-STAT* *hs-TnI* Conversion Guide NPE-5959.1 01/25

