SRČANI TROPONIN I/ (cTnI)

Namjena

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*Heparin pri 90 U/mL dokazano snižava vrijednost cTnI za približno 20 %.

Reference

1. Braunwald, E, *et al.* ACC/AHA 2002 guideline update for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on the Management of Patients with Unstable Angina). 2002. Dostupno na stranici: <http://content.onlinejacc.org/article.aspx?articleid=1130417>.
2. Thygesen K, Alpert JS, White HD, *et al.* Universal definition of myocardial infarction. *Circulation*. 2007.; 116:2634-2653.
3. Anderson JL, Adams CD, Antman EM, *et al.* ACC/AHA 2007 Guidelines for the management of patients with unstable angina/non-ST-elevation myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Writing Committee to Revise the 2002 Guidelines for the Management of Patients With Unstable Angina/Non-ST-Elevation Myocardial Infarction). *Circulation*. 2007.;116:e148-e304.
4. Antman EM, Tanasijevic MJ, Thompson B, *et al.* Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *NEJM* 1996., 335(18): 1342-1349.
5. Galvani M, Ottani F, Ferrini D, *et al.* Prognostic influence of elevated values of cardiac troponin I in patients with unstable angina. *Circulation* 1997., 95: 2053-2059.
6. Morrow DA, Rifai N, Tanasijevic MJ, *et al.* Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes: A thrombolysis in myocardial infarction (TIMI) IIB substudy. *Clin Chem* 2000., 46(4): 453-460.
7. Sabatine MS, Morrow DA, de Lemos JA, *et al.* Multimarker approach to risk stratification in non-ST-elevation acute coronary syndromes: Simultaneous assessment of troponin I, C-reactive protein, and B-type natriuretic peptide. *Circulation* 2002., 105: 1760-1763.
8. Cannon CP, Weintraub WS, Demopoulos LA, *et al.* Comparison of early invasive and conservative strategies in patients with unstable coronary syndromes treated with the glycoprotein IIb/IIIa inhibitor tirofiban. *NEJM* 2001., 344(25): 1879-1887.
9. Morrow DA, Antman EM, Tanasijevic MJ, *et al.* Cardiac troponin I for stratification of early outcomes and the efficacy of enoxaparin in unstable angina: A TIMI-IIB substudy. *JACC* 2000., 36: 1812-1817.
10. Hamm CW, Heeschen C, Goldmann B, *et al.* Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels (CAPTURE Study Investigators). *NEJM* 1999., 340: 1623-1629.
11. Babuin and Jaffe. Troponin: the biomarker of choice for the detection of cardiac injury. *Can. Med. Assoc. J.* 2005.; 173:1191.
12. Hickman *et al.* Cardiac troponin may be released by ischemia alone, without necrosis. *Clin Chem Acta* 2010.; 411: 318-323.
13. Pogledajte "Troponin: What Laboratorians Should know to Manage Elevated Results" na adresi <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/TipsandArticlesonDeviceSafety/ucm109362.htm>.
14. Wu *et al.* NACB Standards of Laboratory Practice: Recommendations for the use of cardiac markers in coronary artery disease. *Clin. Chem.* 1999.; 45:1104.
15. CLSI. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2002.
16. P.J. Cornbleet and N. Gochman, "Incorrect Least-Squares Regression Coefficients in Method-Comparison Analysis," *Clinical Chemistry* 25:3, 432 (1979.).

17. CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline*. CLSI document EP7-A [ISBN 1-56238-480-5]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2002.
18. Clinical and Laboratory Standards Institute (CLSI) *Immunoassay Interference by Endogenous Antibodies; Proposed Guideline*. CLSI document I/LA30-P (ISBN 1-56238-633-6) Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2007.
19. Bjerner et al. Immunometric Assay Interference: Incidence and Prevention. *Clin. Chem.* 2002.; 48:613.
20. Kricka, Interferences in Immunoassays – Still a Threat. *Clin. Chem.* 2000.; 46:1037.
21. Schroff et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res.* 1985.;45:879.
22. Primus et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem.* 1988.; 34:261.
23. Nahm et al. Heteroantibody: phantom of the immunoassay. *Clin. Chem.* 1990.; 36:829.
24. Boscatto et al. Heterophilic antibodies: a problem for all immunoassays. *Clin. Chem.* 1988.; 34:27.

i-STAT is a trademark of the Abbott group of companies.

Abbott Point of Care Inc.
Abbott Park, IL 60064 • SAD

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
Nizozemska

©2018. Abbott Point of Care Inc. Sva prava pridržana. Tiskano u SAD-u