KREATIN KINAZA MB / (CK-MB)

Namjena

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Reference

1. Braunwald, E, et al. ACC/AHA 2002 guideline update for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on the Management of Patients with Unstable Angina). *J Am Coll Cardiol* 2002., 40: 1366-1374.
2. D.W. Moss, A.R. Henderson, "Enzymes" in *Tietz Textbook of Clinical Chemistry – Second Edition*, C.A. Burtis and E.R. Ashwood, eds. (Philadelphia: W.B. Saunders Company, 1994.).
3. Apple FS, Murakami MA. Cardiac Troponin and Creatine Kinase MB Monitoring during In-Hospital Myocardial Reinfarction, *Clin Chem* 2005., 51(2): 460-463.
4. Braunwald E, Antman EM, Beasley JW, Califf RM, Cheitlin MD, Hochman JS, et al. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction. *J Am Coll Cardiol* 2000., 36: 970-1062.
5. Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial infarction defined – a consensus document of the joint European Society of Cardiology/American College of Cardiology Committee for the redefinition of myocardial infarction. *J Am Coll Cardiol* 2000., 36: 959-969.
6. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, et al. It's time for a change to a troponin standard. *Circulation* 2000., 102: 1216-1220.
7. Newby LK, Alpert JS, Ohman EM, Thygesen K, Califf RM. Changing the diagnosis of acute myocardial infarction: implications for practice and clinical investigations. *Am Heart J* 2002., 144: 957-980.
8. Apple FS, Sharkey SW, Falahati A, Murakami MA, Mitha N, Christensen D. Assessment of left ventricular function using serum cardiac troponin I measurements following myocardial infarction. *Clinica Chimica Acta* 1998., 272: 59-67.
9. A.S. Maisel, "Point-of-Care Diagnosis and Management of Myocardial Infarction and Congestive Heart Failure" in *Principles & Practice of Point-of-Care Testing*, G.J. Kost, ed. (Philadelphia: Lippincott Williams & Wilkins, 2002.).
10. Clinical and Laboratory Standards Institute (CLSI). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition*. CLSI document EP9-A2 [ISBN 1-56238-472-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2002.
11. P.J. Cornbleet and N. Gochman, "Incorrect Least-Squares Regression Coefficients in Method-Comparison Analysis," *Clinical Chemistry* 25:3, 432 (1979.).
12. Clinical and Laboratory Standards Institute (CLSI) *Interference Testing in Clinical Chemistry; Approved Guideline*. CLSI document EP7-A [ISBN 1-56238-480-5]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2002.
13. Clinical and Laboratory Standards Institute (CLSI). *Immunoassay Interference by Endogenous Antibodies; Proposed Guidelines*. CLSI document I/LA30-P (ISBN 1-56238-633-6) Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2007.
14. Bjerner et al. Immunometric Assay Interference: Incidence and Prevention. *Clin. Chem.* 2002.; 48:613.
15. Kricka, Interferences in Immunoassays - Still a Threat. *Clin. Chem.* 200.; 46:1037.
16. Schroff et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res.* 1985.; 45:879.
17. Primus et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin. Chem.* 1988.; 34-261.
18. Nahm et al. Heteroantibody: phantom of the immunoassay. *Clin. Chem.* 1990.; 36:829.
19. Boscatto et al. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem.* 1988.; 34:27.

i-STAT is a trademark of Abbott.

Abbott Point of Care Inc.
Abbott Park, IL 60064 • USA

©2026 Abbott Point of Care Inc. All rights reserved. Printed in USA