



# i-STAT/DE

## Customizing Reportable Ranges

### 1 Access the Customization Workspace

#### Notes:

- Reportable Range customization is available on i-STAT/DE versions 2.3 and above.
- Feature is not customizable on the i-STAT 1 analyzer.

### 2 Make sure the "Enable Customization" box has a check mark in it.

Also, make certain that the **Enabled** box is checked for the particular location to which the i-STAT 1 analyzer is assigned.

Location	Enabled
ER	<input checked="" type="checkbox"/>
Lab	<input checked="" type="checkbox"/>

### 3 If the location where this analyzer is assigned has a check mark under the ① Uses Default column, click the alphanumeric code under ② Preferences in the ③ Default Customization Profile column. Otherwise, click the alphanumeric code under ④ Preferences for the specific location to which this analyzer is assigned.

### 4 Once the Preferences window opens, click on the Results tab.

### 5 In the Custom Reportable Ranges section, use the page numbers below the table to find the analyte row which you wish to customize, and click Edit at the right end of that row.

Analyte	Reference Ranges		Action Ranges		Custom Reportable Ranges		
	Low	High	Low	High	Low	High	
Na	138	146	-99999.9	99999.9	-99999.9	99999.9	Edit
K	3.5	4.9	-99999.9	99999.9	-99999.9	99999.9	Edit
Cl	98	109	-99999.9	99999.9	-99999.9	99999.9	Edit

The analyzer software "Default" range values for the analyte will be displayed in the i-STAT/DE Customization Workspace as **-99999.9 for the Low** and **99999.9 for the High**. Use of the default values, allows the analyzer to apply reportable ranges based on the i-STAT Cartridge used for testing.

For reportable range values, see the i-STAT Cartridge Instructions for Use (IFU) and Cartridge and Test Information (CTI) sheet located at [www.pointofcare.abbott](http://www.pointofcare.abbott).

To Customize a reportable range, select the analyte, click ① Edit, enter in the value for the for the **Low** and/or the **High** that is within the reportable range provided in the IFU or CTI sheet. Click on ② Update to save or ③ Cancel to cancel.

**i-STAT ACTk and ACTc Cartridges:** Reportable range is not customizable for these cartridges.


### 6 Click OK and answer OK to the question "Do you want to change Preferences?"

### 7 Dock the analyzer(s) in the downloader which corresponds to the location that the analyzer is assigned. This action will upload the chosen customization features into the analyzer. Repeat step 7 for all analyzers from the same location to be customized. To customize analyzers from other locations for the same features, return to step 2.

# Confirming Custom Reportable Ranges on the i-STAT 1 Analyzer and Display Screens

**1** To view the reportable range on the analyzer, power on the analyzer and press **MENU** once to get to the Administration Menu.


**2** Press **4 -Customization**.




**3** Press **1-View..**




**4** Press **5-Results**.



**5** Press **2-Display Ranges**.



**6** To view displayed Reportable Ranges for each analyte, press the **→** on the analyzer keypad. Blood Gas analytes (pH, PCO<sub>2</sub>, and PO<sub>2</sub>) will display "See IFU" for default and customized Reportable Ranges. See the i-STAT Cartridge IFU or CTI for default Reportable Ranges limits or the i-STAT/DE customization workspace for customized Reportable Range limits.




**NOTE:** Once customized, the reportable range entered will apply to the Patient, Proficiency and Quality Control test pathways on the analyzer but not to the Cal Ver pathway.

## Reportable Range Flags

**Out of Reportable Range - Low**  
A "<" will prefix the result.



**Out of Reportable Range - High**  
A ">" will prefix the result



**Note:** Images are for representational purposes only.

770547-00 Rev.A Revision Date: 24-SEP-2020

i-STAT IS FOR *IN VITRO* DIAGNOSTIC USE ONLY.

For product information and intended use, visit [www.pointofcare.abbott](http://www.pointofcare.abbott)

i-STAT is a trademark of the Abbott. Not all products available in all markets.