



Abbott Point of Care  
400 College Road  
Princeton, NJ 08540

T: +1 609 454 9000  
F: +1 609 419 9370

**DISCONTINUANCE NOTICE**  
**For U.S. Customers using the *i-STAT 1 System***

Product Name	List Number
<b><i>i-STAT EG6+</i> cartridge</b>	03P77-25

**December 2025**

Dear Valued Abbott Point of Care Customer,

This letter contains important information regarding the ***i-STAT EG6+*** cartridge with the *i-STAT 1 System*.

The ***i-STAT EG6+* cartridge** with the *i-STAT 1 System* is intended for use in the *in vitro* quantification of sodium, potassium, hematocrit, pH, partial pressure of oxygen (PO<sub>2</sub>), and partial pressure of carbon dioxide (PCO<sub>2</sub>) in arterial, venous or capillary whole blood.

Abbott Point of Care has made the decision to discontinue the ***i-STAT EG6+*** cartridge for use with the *i-STAT 1 System*. The ***i-STAT EG7+*** and ***CG8+*** cartridges are available on the *i-STAT 1 System* and may be considered as alternative cartridges for your facility. Abbott Point of Care is providing this information in advance to support your planning process and help ensure a smooth transition to an alternative solution that meets your testing needs.




The ***i-STAT EG7+*** and ***i-STAT CG8+*** cartridges have the same intended use as the ***i-STAT EG6+*** cartridge, with the following differences:

- ***i-STAT EG7+*** also includes ionized calcium (iCa) as an additional analyte.
- ***i-STAT CG8+*** also includes ionized calcium (iCa) and glucose as additional analytes.

Abbott Point of Care will stop shipment of the ***i-STAT EG6+* cartridge** on **March 31, 2027**. ***i-STAT EG6+* cartridges** can continue to be used until the specified expiration date on the product label.

## RECOMMENDED ACTIONS:

Please transition to another available *i-STAT* cartridge that offers the range of tests needed on the *i-STAT 1 System*. When transitioning to another *i-STAT* cartridge(s), facilities should follow their own policies and procedures for handling validation testing, if required.

	Discontinued <i>i-STAT</i> cartridge	Available <i>i-STAT</i> cartridges	
Test	<b>EG6+</b> (03P77-25) 	<b>EG7+</b> (03P76-25) 	<b>CG8+</b> (03P88-25) 
Sample Volume	95 µL	95 µL	95 µL
Na	✓	✓	✓
K	✓	✓	✓
iCa		✓	✓
Glu			✓
Hct	✓	✓	✓
Hb*	✓	✓	✓
pH	✓	✓	✓
PCO <sub>2</sub>	✓	✓	✓
PO <sub>2</sub>	✓	✓	✓
TCO <sub>2</sub> *	✓	✓	✓
HCO <sub>3</sub> *	✓	✓	✓
BE*	✓	✓	✓
sO <sub>2</sub> *	✓	✓	✓

\*Indicates a calculated analyte.

## ADDITIONAL INFORMATION:

We appreciate and thank you for your continued support of Abbott and Abbott products. Should you or any of the health care providers you serve have any questions regarding this information, please contact your local support representative or technical services at 1-844-256-9531.

*i-STAT* is for *in vitro* diagnostic use only.

*i-STAT* is a registered trademark of Abbott.

For complete intended use and product information, visit [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott).

11-2025-80 APOC *i-STAT* EG6+ Cartridge on *i-STAT* Stop Shipment Letter US.