



A HIGH-SENSITIVITY TROPONIN ASSAY

# i-STAT hs-TnI CARTRIDGE:

*i-STAT hs-TnI* cartridge is for use with the *i-STAT System* to measure cardiac troponin I (cTnI) in whole blood and plasma to aid in diagnosis of myocardial infarction.

## i-STAT TROPONIN ASSAY CARTRIDGE COMPARISON

RESULTS IN  
≈10 MINUTES



RESULTS IN  
≈15 MINUTES




PRODUCT NAME	i-STAT cTnI	i-STAT hs-TnI
PRODUCT CODE	03P90-25	09P81-25
ANALYTE MEASURED	cTnI	cTnI
SAMPLE SIZE	17 µL	22 µL
SAMPLE TYPE	Venous whole blood or plasma	Venous whole blood or plasma
REPORTABLE RANGE	0.00 – 50.00 ng/mL	2.9 – 1000.0 ng/L
TEST RUN-TIME	Up to 600 seconds (10 minutes)	Approximately 900 seconds (15 minutes)
STORAGE REQUIREMENTS	Refrigerated at 2-8°C (35-46°F) until expiration date. Room temperature at 18-30°C (64-86°F) for 2 weeks	Refrigerated at 2-8°C (35-46°F) until expiration date. Room temperature at 18-30°C (64-86°F) for 2 weeks
CARTRIDGE BOX DIMENSIONS	3.75" (W) x 2" (H) x 9.625" (L) 25 cartridges per box	3.75" (W) x 2" (H) x 9.625" (L) 25 cartridges per box
CARTRIDGE POUCH/PACK DIMENSIONS	3.7" (W) x 1.95" (H) x 0.37" (D)	3.7" (W) x 1.95" (H) x 0.37" (D)
CONTROLS	Liquid: 06P17-09 Level 1   06P17-10 Level 2   06P17-11 Level 3	Frozen: 06P17-21 Level 1   06P17-22 Level 2   06P17-23 Level 3

i-STAT


# i-STAT TROPONIN ASSAY CARTON PACKAGING ARTWORK COMPARISON

03P90-25

09P81-25





TOP  
PANEL




FRONT  
PANEL




BACK  
PANEL



LEFT & RIGHT  
SIDE PANELS



BOTTOM  
PANEL



The *i-STAT* *hs-TnI* cartridge with the *i-STAT* System is intended to be used as an aid in the diagnosis of myocardial infarction (MI).  
For *in vitro* diagnostic use only. This material is intended for use in Canada.  
For intended use and complete product information, visit [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott)

©2025 Abbott. All rights reserved. *i-STAT* is a trademark of Abbott.  
Abbott | Point of Care Diagnostics | 400 College Road East | Princeton, NJ 08540 USA | 609 454 9000 | [www.globalpointofcare.abbott](http://www.globalpointofcare.abbott) | *i-STAT* *hs-TnI* vs *cTnI* Product Differences Flashcard | CA | 25002454.1 05/25

