



i-STAT G3+ CARTRIDGE (WHITE): NOW FDA CLEARED

i-STAT G3+ cartridge is now FDA cleared for use with the *i-STAT 1 analyzer* as part of the *i-STAT 1 System* to measure pH, partial pressure of oxygen (PO₂), and partial pressure of carbon dioxide (PCO₂) in patients of all age ranges, neonates through adults.



DISCONTINUED




NOW FDA CLEARED

CARTRIDGE TITLE	<i>i-STAT G3+</i> CARTRIDGE (BLUE)	<i>i-STAT G3+</i> CARTRIDGE (WHITE)
PRODUCT CODE	03P78-50**	03P78-26
ASSAYS MEASURED	pH, PO ₂ , PCO ₂	pH, PO ₂ , PCO ₂
SAMPLE SIZE	95 uL	95 uL
SAMPLE TYPE	ARTERIAL AND VENOUS	ARTERIAL, VENOUS AND CAPILLARY
REPORTABLE RANGE	pH 6.5-8.2 PO ₂ 5-800 mmHg PCO ₂ 5-130 mmHg	pH 6.5-7.8 PO₂ 5-700 mmHg PCO ₂ 5-130 mmHg
TEST RUN TIME	~130-200 SECONDS	~130-200 SECONDS
STORAGE REQUIREMENTS	REFRIGERATED AT 2-8°C (35-46°F) UNTIL EXPIRATION DATE. ROOM TEMPERATURE AT 18-30°C (64-86°F) FOR 2 MONTHS	REFRIGERATED AT 2-8°C (35-46°F) UNTIL EXPIRATION DATE. ROOM TEMPERATURE AT 18-30°C (64-86°F) FOR 2 MONTHS
CARTIRIDGE BOX DIMENSIONS	2.75" (H) x 8.00" (W) x 3.81" (D)	2.75" (H) x 8.00" (W) x 3.81" (D)
CARTRIDGE POUCH DIMENSIONS	2.64" (H) x 3.74" (W)	2.64" (H) x 3.74" (W)


***i-STAT G3+* Product Code 03P78-50 was discontinued on August 31, 2023

i-STAT G3+ PACKAGING ARTWORK DIFFERENCES

i-STAT G3+ 03P78-50



i-STAT G3+ 03P78-26



For *in vitro* diagnostic use only. Not all products available in all regions.

This material is intended for a U.S. audience only.

©2023 Abbott. All rights reserved. *i-STAT* is a trademark of Abbott.

Abbott | Point of Care Diagnostics | 400 College Road East | Princeton, NJ 08540

609 454 9000 | www.globalpointofcare.abbott

i-STAT G3+ Product Differences Flashcard | 4966.REV1 09/13

